# $\alpha_1$ -ac. glycoprotein at

*Turbidimetric method*

<b>REF CT31350</b> 2 x 40 mL <b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>REF CT31352</b> 6 x 40 mL <b>CONTENTS</b> R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

$\alpha_1$ -acid glycoprotein at is a quantitative turbidimetric assay<sup>1</sup> for the measurement of  $\alpha_1$ -acid glycoprotein ( $\alpha_1$ -AG) in human serum or plasma. Anti-human  $\alpha_1$ -AG antibodies form insoluble complexes when mixed with samples containing  $\alpha_1$ -AG. The light scattering of the immunocomplexes depends upon the  $\alpha_1$ -AG concentration of the patient sample, and can be quantified by comparison from a calibrator of known  $\alpha_1$ -AG concentration.

## REAGENTS

**R1  $\alpha_1$ -acid glycoprotein at.** Goat antibodies anti-human  $\alpha_1$ -acid glycoprotein, tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.95 g/L.

**Precautions:** The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

## PREPARATION

**R1.** Ready to use.

### Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each  $\alpha_1$ -acid glycoprotein Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

## STORAGE AND STABILITY

- Store at 2-8°C.

The reagent is stable until the expiration date stated on the label, when stored **tightly closed** at 2-8°C and contaminations are prevented during their use. Do not use the reagent after the expiration date.

- Presence of particles, turbidity and the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

## SAMPLE COLLECTION

Fresh serum or, EDTA or heparinized plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

$\alpha_1$ -AG in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

## INTERFERENCES

- Bilirubin (10 mg/dL), hemoglobin (2 g/L), rheumatoid factors (200 IU/mL), do not interfere.
- Lipemia ( $\geq 1.25$  g/L), may affect the results. Other substances may interfere<sup>6,7</sup>.

## ADDITIONAL EQUIPMENT

- Analyzer LIDA.

## AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.

A reagent blank should be run daily before sample analysis.

## REFERENCE VALUES

Adults<sup>2</sup>: 50 – 120 mg/dL.

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS<sup>3-5</sup>

$\alpha_1$ - acid glycoprotein (orosomucoide) is an acute phase protein<sup>3</sup> shown to exert a regulatory, dampening influence on the inflammatory cascade, thereby protecting against tissue damage from excessive inflammation.

$\alpha_1$ - acid glycoprotein is synthesized by hepatic parenchymal cells, but granulocytes and monocytes may also contribute significantly to plasma levels in sepsis. It is classified as a lipocalin, a group of proteins that binds lipophilic substances as progesterone and related hormones<sup>5</sup>.

Plasma concentration increases 3 to 4-fold in most conditions associated with inflammation or tissue necrosis, and may be one of the most reliable indicators of clinical activity of ulcerative colitis<sup>4</sup>. Levels also are increased by glucocorticoid effect, either endogenous (Cushing's syndrome) or exogenous, along with heptaglobin and prealbumin levels. Synthesis and plasma levels are decreased by estrogens<sup>5</sup>.

Patients with proteinuria,  $\alpha_1$ -acid glycoprotein is preferentially excreted in the urine.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## BIBLIOGRAPHY

- Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
- Dati F et al. *Eur J Clin Chem Clin Biochem* 34:517-520 (1996).
- Kushner, Mackiewicz A. *Disease Markers* 1-11 (1987).
- Birger K et al. *Scand J Gastroent* 11: 177-183 (1976).
- Tietz *Textbook of Clinical Chemistry*, 3<sup>rd</sup> Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
- Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
- Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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## REFERENCIAS

1. Price CP et al. Ann Clin Biochem 20: 1-14 (1983).
2. Dati F et al. Eur J Clin Chem Clin Biochem 34:517-520 (1996).
3. Kushner, Mackiewicz A. Disease Markers 1-11 (1987).
4. Birger K et al. Scand J Gastroent 11: 177-183 (1976).
5. Tietz Textbook of Clinical Chemistry, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
6. Young DS. Effects of drugs on clinical laboratory tests. 3th ed. AACC Press (1997).
7. Friedman and Young. Effects of the disease on clinical laboratory tests, 3th ed. AACC Press, 1997.