

LINEAR

Your Reliable Partner in IVD

since
1994

Gastroenterology


Rapid Results • Open System • Easy to use • High Throughput

1 sample
muestra **3** tests

iFOB • iTRA • iCALPRO




Our Stool Collection tube



- The buffer inside keeps stable the sample during 7 days at room temperature (20°-25°) and 10 days at 2-8°c
- Barcoded
- Open System: applicable in almost all the Automated Analyzers and in our Kroma iT series

Sample Collection




Our dedicated instruments

KROMA iT



- Throughput: 80 tests per hour
- 27 sample positions
- 9 reagent positions
- Washing station
- Barcode reader
- LIS bi-directional connection
- Calibration graphics
- Cap piercing for SCT tubes

KROMA PLUS iT



- Throughput: 150 tests per hour
- 49 Sample positions
- 35 Reagent positions
- Washing station
- Barcode reader
- LIS bi-directional connection
- Calibration graphics
- Cap piercing for SCT tubes